

# BITREX<sup>®</sup> Fit Test Kit 0102



## FIT TEST RESPIRATOR WEARERS

Test the face to facepiece seal of Moldex<sup>®</sup> particulate respirators. Check that respirators are being worn correctly and that they fit satisfactorily. Ensure that users are properly protected.

## EASY TO USE

The Moldex<sup>®</sup> Fit Test Kit contains everything needed to conduct a qualitative fit test. Written instructions provide easy-to-use, step-by-step instructions. It's a simple pass/fail test. Taste sensitivity to BITREX<sup>®</sup> is checked. The wearer dons a Moldex<sup>®</sup> particulate respirator, the fit test hood is put on and the fit test solution is sprayed into the hood.

## APPROVED FIT TEST METHOD

OSHA requires that a user be fit tested before wearing a respirator. The Moldex<sup>®</sup> BITREX<sup>®</sup> Fit Test Kit can help you meet respirator program fit test requirements established by OSHA 1910.134.

## MOLDEX<sup>®</sup> TECHNICAL SUPPORT

Moldex<sup>®</sup> provides technical assistance including help with fit testing, a Train-the-Trainer Program, a manual on How To Set Up Your Respiratory Protection Program, plus a Toll-Free Technical Help line. Just call 800/421-0668 ext. 512/550 or 310/837-6500 ext. 512/550.

## FEATURES

- Includes everything necessary to conduct a BITREX<sup>®</sup> Qualitative Fit Test.
- Easy-to-follow, step-by-step instructions.
- Helps satisfy OSHA respirator program requirements.





## 0102 BITREX® QUALITATIVE FIT TEST KIT

0201 Fit Test Hood (1)  
0301 Nebulizer #1 (Threshold Screening for Taste Sensitivity)  
0302 Nebulizer #2 (Fit Test)  
0503 Sensitivity Screening Solution (2.5 mL per Ampule)  
0504 Fit Test Solution (2.5 mL per Ampule)  
0702 Qualitative Fit Test Instructions

### NOTES:

1. Use with any Moldex® N, R or P class disposable respirator. This fit test kit may also be used with any Moldex® half mask respirator by using any N, R or P class filter.
2. You should get approximately 3 – 5 fit tests per applicator tube (fit test cycle) or 18 – 30 per supply of BITREX® solution in this kit. (Reorder #0503 & #0504)

### WARNING TO USER

1. Before use, the user must first be trained by the employer in proper respirator use, in accordance with applicable safety and health standards, for the contaminant and exposure level in the assigned work area.
2. The particulates and other contaminants which can be dangerous to your health include those which you cannot see, taste, or smell.

### LIMITED WARRANTY IMPORTANT NOTICE TO PURCHASER

This limited warranty is made in lieu of the warranties of merchantability, fitness for particular purposes and all other warranties, express or implied. There are no other warranties which extend beyond the description on the face hereof. The physical standards and specifications of Moldex will be met by products sold. **Exclusive Remedies:** damages for the breach of this limited warranty are limited to the replacement of such quantity of Moldex products proved to be defectively manufactured. Except as provided above, Moldex shall not be liable or responsible for any loss, damages, or liability, direct, indirect, incidental, special or consequential, arising out of a sale, use or misuse, or the inability to use products by the user.

### DISTRIBUTED BY:

#### MOLDEX-METRIC, INC.

10111 W. Jefferson Blvd. Culver City, CA 90232  
TEL: 800/421-0668 or 310/837-6500 FAX 310/837-9563  
www.moldex.com e-mail: sales@moldex.com

#### CANADIAN CUSTOMER SERVICE

TEL: 310-837-6500, Ext. 517  
FAX: 310-837-9563

Moldex, Ideas that wear well, Dura-Mesh and Ventex are registered trademarks.  
The mesh design is a trademark of Moldex-Metric, Inc.  
U.S. Patents #D285,374 and #4,850,374.  
Made in the U.S.A. 9700-790 REV F 01/09



Printed on Recycled Paper



**MOLDEX**  
Ideas that wear well.®